



Timing function



Suitable for a variety of containers



Adjustable speed




3D mixing

CSK - 3D



3D mixing - Adjustable speed



10° Inclination meets a wide range of applications



Non-slip silicone tray surface



Suitable for a variety of containers



TECHNICAL DATA

Model	CSK - 3D
Type	3D Oscillating
Display	LCD
Tit Angle	10°
Max Load	0.8 kg
Speed	10 - 80 rpm
Timing	1min - 8h
Allowed Ambient Humidity Temp.	5 ~ 40 °C 80%RH
Voltage	DC 12V
Power	10W
Dimension (mm)	184*184*130