



Timing function



Stepless speed regulation, wide speed range




LCD Display



Durable brushless motor

CMMS4 Pro / CMMS6 Pro / CMMS9 Pro



TECHNICAL DATA

| Model | CMMS4 Pro | CMMS6 Pro | CMMS9 Pro |
|---------------------------|----------------------|---------------|---------------|
| Control | Asynchronous Control | | |
| Number of Stirring Points | 4 | 6 | 9 |
| Max Capacity | 4*500ml | 6*500ml | 9*500ml |
| Working Plate Size | 4-(110*110mm) | 6-(110*110mm) | 9-(110*110mm) |
| Motor | Brushless motor | | |
| Speed | 50-1500rpm | | |
| Bar Size (mm) | Z30 | | |
| Timing | 1min-19h 59min | | |
| Voltage | AC100-240V | | |