



Memory program



Multi rotors optional



Temperature protection & overshock protection



Double electronic lock

LC500-8 / LC500-12 / LC500-6 / LC500-18 / LC500-24



TECHNICAL DATA

Model	LC500-8	LC500-12	LC500-6	LC500-18	LC500-24
Rotor Capacity	Centrifuge tube: 15 ml*8	Centrifuge tube: 15ml*12	Centrifuge tube: 15 / 50ml*6	Blood tube: 2 / 5 / 10ml*18	Blood tube: 2 / 5 / 10ml*24
	Blood tube: 2 / 5 / 10ml*8	Blood tube: 2 / 5 / 10ml*12	Blood tube: 2 / 5 / 10ml*6		
Speed	500 - 4000 rpm				
Max RCF	1967 x g	2075 x g	2146 x g	2182 x g	2164 x g
Semidiameter of Rotor	11 cm	11.6 cm	12 cm	11.2 cm	12.1 cm
Timing Range	30s - 59min 50s				
Protection Function	Lid opening protection, overtemperature protection, overshock protection				
Power	60W	60W	80W	80W	80W
Power Supply	AC 110V - 220V, 50/60 Hz				