



LAB CENTRIFUGE Digital Mini Centrifuge

CYTON
LABORATORY INSTRUMENTS

B

www.cytosci.com



LCD Display



Multi-rotor, big working capacity



Run smoothly, precise separation



Speed adjustable and timing

MC - 5Pro / MC - 7 Pro

Power Switch



Display




Safety Lock


Rotor Elastic Buckle



Rotor



Rotor 0.2/0.5/1.5/2ml*12



Rotor 40ml capillary*12

TECHNICAL DATA

Model	MC - 5 Pro	MC - 7 Pro
Control Mode	Digital Display & Button	
Display	LCD Display	
Max Relative Centrifugal Force	1680 x g	3286 x g
Max Rotation Speed	5000 rpm $\pm$ 5%	7000 rpm $\pm$ 5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 Centrifuge tube 0.5/1.0ml*12 Micro blood tube	
Timing	10sec - 99min	
Power	DC 24V	