



Multiple protections, more secure



Two calibration modes



High-efficiency halogen heating lamp



Multi-device printing, standard RS232 interface

MBA1102 / MBA1102a / MBA1103 / MBA1103 (HBM)



TECHNICAL DATA

Model	MBA1102	MBA1102a	MBA1103	MBA1103(HBM)
Max Capacity	110g			
Weighing Precision	0.01 g	0.005 g	0.01 g	0.001 g
Sensor	Load Cell	Load Cell	Load Cell	HBM
Temp.Sensor	PT-100-2			
Temp.Range	40 °C-160 °C			
Temp.Setting / Display Accuracy	1 °C			
Pan Size	φ 96 mm			
Moisture / Testing Range	90.00% - 10.00%	95.00% - 5.00%	100.00% - 0.00%	100.00% - 0.00%
Moisture / DRY% Readability	0.1%	0.05%	0.01%	0.01%
Operate. Temp	5 °C - 35 °C			
Language	CN / EN			