



Stepless speed regulation  
Wide speed range



Press to start / Compact body



High quality and colorful silicone



Brushless motor / Durable

VM - 300 / VM - 300S

Oscillator

Speed Control Knob

Small size, easy to grasp




300 rpm fast mixing



Can be used for shaking Test tubes,  
Centrifuge tube, Eppendorf tubes,  
with significant mixing effect.

Start protection structure



Technical innovation, special starting protection device, longer service life

Zinc alloy base



The weight is about three times that of aluminum alloy, and the body is more stable.

## TECHNICAL DATA

Model	VM-300	VM-300S
Control Mode	Constant speed	Knob
Motor	Brushless motor	
Amplitude	4.5mm (Circular oscillation)	
Maximum Sample Processing Capacity	50ml	
Speed Range	3000rpm	0-3000rpm

Shell Protection Level	IP42	
Power	12W	
Voltage	100-240V, 50 / 60Hz	
Product Size	93*93*72mm	100*93*72mm
Allowable Ambient Temperature	5-40°C, 80%RH	
Gross / Net weight	0.55 / 0.36kg	0.56 / 0.37kg