



# LAB CENTRIFUGE Mini Centrifuge - (K) Series

CYTON  
LABORATORY INSTRUMENTS

B

www.cytonsci.com



Stop running after opening lid



Low noise

MC-4K / MC-7K / MC-10K / MC-12K

Visible transparent top lid

Rotor

MC-4K

Power Switch

The silicon feet

Rotor 0.2ml\*16 + 1.5 / 2.0ml\*8

Rotor 1.5 / 2.0ml\*8

## TECHNICAL DATA

| Model           | MC-4K                    | MC-7K       | MC-10K        | MC-12K      |
|-----------------|--------------------------|-------------|---------------|-------------|
| Max Speed       | 4000rpm                  | 7000rpm     | 10000rpm      | 12000rpm    |
| Max RCF         | 984xg                    | 3013xg      | 6149xg        | 8855xg      |
| Rotor Capacity  | 0.2ml*16 + 1.5 / 2.0ml*8 |             | 1.5 / 2.0ml*8 |             |
| Power           | 15W                      | 20W         | 30W           | 40W         |
| Dimensions (mm) | 160*135*110              | 160*135*110 | 160*135*110   | 160*135*110 |
| Weight (kg)     | 0.8                      | 0.8         | 0.8           | 0.8         |