



Timing function



Large handling capacity


VM - 600max



Operation Panel

Time Control Knob

Speed Control Knob



Set time

Set speed

Rotate: Set Time / Tap: Start - Pause


Rotate: Set the RPM / Tap - Start - Pause

Open type fixture




Open type modules, more convenient to take samples greatly improve the operation convenience of the experiment.

Small body - large handling capacity



Compact and beautiful products, large capacity, the Max processable sample: 15*16+50*1*25

Multiple modules optional



Multiple modules are available simple handling of various specifications. Experiments are more convenient and quick.

Widely use



Widely used in medical experiment, life sciences, chemical pharmaceuticals, agricultural scientific research, etc.

TECHNICAL DATA

Model	VM - 600Max
Motor Type	Brushless motor
Horizontal Gyration Amplitude	4mm
Speed	1500 - 3000 rpm
RPM Increment	10rpm
Timing	1min - 8h
Handling Range	Model (VM601) Centrifuge tube: 50ml*25+15ml*16 Model (VM602) Centrifuge tube: 15ml*49 / Blood tube: 10ml*49 Model (VM603) Blood tube: 2/5ml*49 Model (VM604) Centrifuge tube: 1.5/2ml*81
Allowed Ambient Humidity Temp.	5 ~ 40 °C 80%RH
Voltage	AC 100 - 220V, 50/60Hz
Power	30W
Dimension (mm)	355*260*145