



Ultra-low body, stable operation



Brushless motor with strong power



Various fixtures, large handling capacity



Touch / Run, easy to switch

VM - 500S / VM - 500Pro




Suitable for mixed operation of microplate




The unique plate is suitable for 96-well microplate, mixing is more convenient.

Multiple adaptive modules



One machine can be equipped with 8 modules, suitable for single tube, multi-tube, flask, etc., could meet your different mixing requirements.



TECHNICAL DATA

Model	VM-500S	VM-500Pro
Operating Mode	Touch / Run	
Control Mode	Knob	Digital display & Button
Display	No	LCD
Motor	Brushless motor	
Amplitude	4.5 mm (Circular oscillation)	

Model	VM-500S	VM-500Pro
Speed Range	Touch: 0-3000 rpm Run: 200-3000 rpm	Touch: 0-3000 rpm Run: 100-3000 rpm
Timing / Time	No	1min - 8h
Shell Protection Level	IP21	
Voltage	20W	
Power	DC12V	