



Sterilization only takes 5-7 seconds



Can be used in an anaerobic environment



No open-flame, more safe



Easy to clean

DS-900 / DS-900S

Heating Zone

Heating hole

DS900

DS900S

Power Switch

Power Light



TECHNICAL DATA

Model	DS - 900	DS - 900S
Power	150W	1240W
Max. Temperature Inside the Heating Hole	825 °C ± 50 °C	
Heating Temperature	Two temperatures available 400 °C / 850 °C	
Max. Diameter of Sterilizer Goods (mm)	Φ14	Φ35
Length of Heating Zone (mm)	150	200
Warm-up Time	10 min	