



LCD Display



Electronic lock



Emergency Switches



Memory program

LC400 / LC400 - PRO



TECHNICAL DATA

Model	LC400	LC400 Pro
Speed Range	500 - 4000 rpm	500 - 5000 rpm
Max Relative Centrifugal Force	1717 x g	2683 x g
Rotor Capacity	15 ml*4 Centrifuge tube 2 / 5 / 10 ml*8 blood tube	
Increment	100 rpm	
Timing Range	30s - 59min 50s	
Noise	≤60dB (A)	
Voltage	AC100 - 240V	
Power	35W	45W
Dimensions (mm)	265*260*190	